

Bigelow | Laboratory for
Ocean Sciences

Single cell genomics workflow

Greg Gavelis, Julia Brown, Ramunas Stepanauskas

Bioinformatics of Microbial Single Cells
Tuesday April 5th, 2022




Microbial single cell genomics workflow




As Ramunas explained this morning, the SCGC process begins with sample collection and lysis, and results in terabytes of data


Microbial single cell genomics workflow




As Ramunas explained this morning, the SCGC process begins with sample collection and lysis, and results in terabytes of data




As Ramunas explained this morning, the SCGC process begins with sample collection and lysis, and results in terabytes of data




As Ramunas explained this morning, the SCGC process begins with sample collection and lysis, and results in terabytes of data



As Ramunas explained this morning, the SCGC process begins with sample collection and lysis, and results in terabytes of data



As Ramunas explained this morning, the SCGC process begins with sample collection and lysis, and results in terabytes of data



Many bioinformatic steps are involved in cleaning, assembling, and annotating the genomes. This presentation will be a brief overview of the different steps involved. When space allows, I include the actual commands that we use at each step of the pipeline.


How we generate the final products

& where to find them

Our BATS datasets has 377 samples, with names like:

AH-141-C18
AH-141-D10
AH-141-I04

etc



Each sample gets its own namesake folder. That folder contains the final genome assembly (e.g. AH-141-C18_contigs.fasta), and also three directories, where the products of varies pipeline steps are deposited. (Note that we don't include all of the intermediates from read processing, since fastq files are very large.)


1. Report Read Quality

PIPELINE

FastQC



```
fastqc -q {forward.fastq} {reverse.fastq}
```



Before processing, we review the quality of the raw reads, using the program `fastqc` to output html report files, like above. These are included in the QC subfolder. Note that the quality of the basecalls tends to decline as we reach the end of the 150 bp read.

2. Trim Reads

PIPELINE

FastQC



Trim Reads

Remove low quality 5' and 3' bases.

```
trimmomatic PE -phred33 \  
    {Forward.fastq} \  
    {Reverse.fastq} \  
    LEADING:0 TRAILING:5 SLIDINGWINDOW:4:15 MINLEN:36
```

For our next step, we use the program trimmomatic to remove these low quality reads. Reads under 36 bp are discarded.

3. Remove low-complexity reads

- Custom script in python

PIPELINE

FastQC



Trim Reads

Filter Low-complexity



Step 3, we remove low-complexity reads include (e.g. homopolymer runs "AAAAAA" and repeats "ATATATA"). If not removed, low complexity sequences can cause misassemblies


4. Normalize reads

For computational efficiency, downsample readpairs with over-represented kmers.

- (our kmers are 21bp subsequences.)

```
kmernorm -k 21 -t 30 -c 3 \  
{readpairs.fastq} > {ID}_normalized_pe.fastq
```

PIPELINE



Kmernorm is much faster than the the read normalization program that comes built into the assembly software. Normalization allows assembly to run ~10X faster and with no significant decline in assembly quality.

5. Remove contaminant reads

- Custom python script
- BWA (Burroughs-Wheeler Aligner)
- Align to known contaminant genomes

PIPELINE

FastQC



Trim Reads

Filter Low-complexity

Normalize




Decontaminate



For the first of three decontamination steps, we remove reads that align to nucleotides in a custom database of genomic contaminants. We obtained these contaminants by sequencing 'empty' (i.e. negative control) plates in the clean lab.

6. Assemble Genome


- **SPAdes** uses de-bruijn graphs (across various kmer sizes) to assemble contigs.
- Unlike other assemblers, it *doesn't rely on read-coverage* (We expect uneven coverage due to MDA)



SPAdes is unparalleled for the “de-novo” (i.e. reference-free) assembly of DNA obtained from single cells. Critically, it assumes that coverage of input reads will be highly uneven (due to amplification bias introduced during Multiple Displacement Amplification).

6. Assemble Genome

```
spades.py -o {ID}_all_contigs.fasta \  
    --careful \  
    --sc \  
    --phred-offset 33 \  
    {reads.fastq}
```




We make sure to run spades in single-cell mode, so that it does not take read coverage into account. This output assembly is considered raw and likely has some areas that need to be processed. The next several slides deal with matters of QC, that have been carefully optimized.

7. Trim contigs

- Most misassemblies occur at ends of contigs (200 bp)

8. Length filtering

- Discard assemblies under 2000 bp.





9. Filter contaminant contigs

- Via **blastn** against a db of known contaminants.

Contaminant contigs are in:

{ID}_contaminated_contigs.fasta




Again, sample nucleotides are aligned to known contaminants (this time with **blastn**.)
All these potential contaminants are retained for the user to review, if desired.

9. Filter contaminant contigs

- Via **blastn** against a db of known contaminants.



Contaminant contigs are in:

{ID}_contaminated_contigs.fasta



- The final cleaned assembly is in:

{ID}_contigs.fasta



10. Preliminary SSU-based classification

- **Blastn** contigs against reference database (Silva) to detect putative 16s.

-> {ID}_SSU.fasta




- **Samtools** aligns 16s hits to closest reference sequences.
Custom python script assigns Linnean classification.

-> {ID}_SSU.tsv

Example:

k__Bacteria;p__Proteobacteria;c__Gammaproteobacteria;o__E01-9C-26_marine_group;f__?;g__?;s__?




This is only one classification approach. In coming lectures we'll talk about how multigene approaches are also used.

12. Genome annotation with Prokka

- **Prokka** bundles a number of programs to annotate genomic features.

- CDS: Coding sequences (Prodigal)
- rRNA (RNAmmer)
- tRNA (Aragorn)
- non-coding RNA (Infernal)
- Signal leader peptides (SignalP)




In coming lectures, we'll also discuss other annotation programs, like DRAM.

12. Prokka outputs



- Prokka Output files:

- **.ffn** – Fasta of all genomic features
- **.faa** – Proteins, i.e. translated CDS
- **.gbk/.gff** – Files of seqs + annotations
- **./tsv.tbl** – Feature tables
- **.txt** – Counts of each feature type



12. Where are prokka outputs stored?

- 'Prokka' dir – Prokka run on prokadb
- 'Swissprot' dir – Prokka run on prokadb + swissprot



Lastly, there is a **stats sheet** to summarize the metadata of the entire sequencing plate....



You'll take a look at that sheet this afternoon, during the interactive session.

Future Improvements:

- GTDB-Tk for better taxonomic classification
- DRAM for deeper metabolic annotation
- And we're always open to suggestions!

